

BAUSCH+LOMB

ANTERIOR

Product Catalog



Stellaris[®] Elite
Vision Enhancement System

Table of Contents

Stellaris Elite®	2
Phaco Packs	4
Phaco Needles and Sleeves	12
Irrigation/Aspiration	18
Instruments	28
Accessories	34



Stellaris Elite®

Stellaris Elite® Anterior Deluxe System
Comes with Stellaris Elite® foot pedal



Stellaris Elite® Foot Control Laser-Integrated



Stellaris® Attune Ultrasonic Handpiece
BL3170
with SureLock



Needle and sleeve not included.

Stellaris® I/A Handpiece
BL3160
with SureLock




Needle and sleeve not included.

I/A Handpiece
DP9721
Titanium



DigiFlow™ Tubing
D4600A
Air exchange line



Digital Media System (DMS)
BL2300HD
Video overlay system



Phaco Packs

Adding convenience and efficiency to procedures, our Phaco packs are configured for your procedural needs. Find the perfect fit for your practice and choose from the selection of components that surgeons rely upon for their quality and dependability.

Easily reordered to keep up to date with your workflow, Phaco packs are a convenient way to add efficiency to your practice. The sterilized, single-use components display the high quality engineering and craftsmanship that surgeons have come to expect from Bausch + Lomb Surgical.



Basic Vacuum Gravity Packs

Stellaris® Basic Vacuum Phaco Pack

BL5111
A 1 DP5501 light blue sleeve



StableChamber® Vacuum Gravity Packs

Stellaris® StableChamber® Vacuum Pack

BL5112
Pack Includes:
A 2 DP5501 light blue sleeves




Stellaris® StableChamber® Vacuum Pack with MICST™ Needle

BL5113
Pack Includes:
A 1 BL3318S purple straight needle
B 2 BL3118 MICST™ sleeves



Stellaris® StableChamber® Vacuum Pack with MICST™ Angled Needle


BL5113A
Pack Includes:
A 1 BL3318AS purple angled needle
B 2 BL3118 MICST™ sleeves



Pack come with the items below.
1 Vacuum cassette and tubing
2 Test chamber
3 Needle wrench




Packs come with the items below.
1 Vacuum cassette and StableChamber® tubing
2 Screen drape
3 Tray cover drape
4 Remote control drape
5 Test chamber
6 Needle wrench
7 Plug



Premium Vacuum Gravity Packs


Stellaris® Premium Vacuum Phaco Pack
BL5110
Pack Includes:
A 2 DP5501 light blue sleeves



Stellaris® Premium Vacuum Phaco Pack with MICST™ Straight Needle
BL5114
Pack Includes:
A 1 BL3318S purple needle
B 2 BL3118 MICST™ sleeves




Stellaris® Premium Vacuum Phaco Pack with MICST™ Angled Needle
BL5114A
Pack Includes:
A 1 BL3318AS purple needle
B 2 BL3118 MICST™ sleeves



Premium Vacuum Gravity Packs


2.6 to 2.8mm Vacuum Phaco Pack
BLV2260000
Pack Includes:
A 2 DP5511 dark blue sleeves



2.2 to 2.4mm Vacuum Phaco Pack
BLV3220000
Pack Includes:
A 2 DP5531 green sleeves




1.8 to 2.0mm Vacuum Pack
BLV4180000
Pack Includes:
A 2 BL3118 MICST™ sleeves



Packs come with the items below.
1 Vacuum cassette and tubing
2 Screen drape
3 Tray cover drape
4 Remote control drape
5 Test chamber
6 Needle wrench
7 Plug




Packs come with the items below.
1 Vacuum cassette and tubing
2 Screen drape
3 Tray cover drape
4 Test chamber
5 Needle wrench
6 Plug



Basic Adaptive Fluidics™ Packs


Adaptive Fluidics™ Basic Pack with MICST™ Sleeves

- BL5115-4
- Pack Includes:**
- A 2 BL3118 MICST™ sleeves
 - B Adaptive Fluidics cassette and tubing



Adaptive Fluidics™ Standard Pack with MICST™ Sleeves and Needle

- BL5115-2
- Pack Includes:**
- A BL3420S MICST™ blue straight needle
 - B 2 BL3124S yellow sleeves
 - C Adaptive Fluidics cassette and tubing with Pulse Reflux




- Above packs also include:**
- 1 Screen drape
 - 2 Tray cover drape
 - 3 Remote control drape
 - 4 Test chamber
 - 5 Bottle spike
 - 6 Needle wrench
 - 7 Plug



Premium Adaptive Fluidics™ Packs


Adaptive Fluidics™ Premium Pack

- BL5115-1
- Pack Includes:**
- A BL3420S MICST™ blue straight needle
 - B 2 BL3124S yellow sleeves



Adaptive Fluidics™ Premium MICST™ Pack

- BL5115-3
- Pack Includes:**
- A BL3420S MICST™ blue straight needle
 - B 2 BL3118 MICST™ sleeves



- Packs come with the items below.**
- 1 Adaptive Fluidics™ cassette and tubing with StableChamber® and Pulse Reflux
 - 2 Screen drape
 - 3 Tray cover drape
 - 4 Remote control drape
 - 5 Test chamber
 - 6 Needle wrench
 - 7 Bottle spike
 - 8 Plug



Phaco Needles and Sleeves

Designed for ultimate precision, manufactured from high grade materials, and available in a wide range of sizes, our selection of Phaco needles and sleeves can suit a variety of needs for many types of practices.

Choose from options including a range of angles, straight needles, MircoFlow®, MicroFlow® Plus, and Cortex with thin tips and various colored threads also available for easy identification in the operating room.



Phaco Needles Single-Use, Sterile 6/box

30° MICST™ 2.0 Angled Needle
BL3420AS
RIS: 2.0mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



45° 19G Straight Needle
DP8145S
RIS: 3.00-3.20mm OD: 1.70mm
ID: 0.9mm



30° MicroFlow® 2.2 Straight Needle
DP8730S
RIS: 2.2mm OD:0.90mm
ID:0.70-0.50mm



30° MICST™ 2.0 Straight Needle
BL3420S
RIS: 2.0mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



30° MicroFlow® Straight Needle
DP8230S
RIS: 2.50-2.65mm OD: 1.07mm
ID: 0.90-0.51mm



30° Vortex 20-21G Straight Needle
DP8930S
RIS: 2.40mm OD:0.84mm
ID:0.60mm



30° MICST™ Angled Needle
BL3318AS
RIS: 1.8-2.00mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



45° MicroFlow® Straight Needle
DP8245S
RIS: 2.50-2.65mm OD: 1.07mm
ID: 0.90-0.51mm



30° Thin Tip 20G Angled Needle
DP8430AS
RIS: 2.65-2.75mm OD: 0.90mm
ID: 0.67mm



30° MICST™ Straight Needle
BL3318S
RIS: 1.8-2.00mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



30° MicroFlow® Plus Straight Needle
DP8330S
RIS: 2.75-2.85mm OD:1.17mm
ID:0.98-0.61mm



30° Thin Tip 20G Straight Needle
DP8430S
RIS: 2.65-2.75mm OD: 0.90mm
ID: 0.67mm



30° 19G Straight Needle
DP8130S
RIS: 3.00-3.20mm OD: 1.70mm
ID: 0.90mm



Phaco Needles - Non-Sterile, Limited Use 6/box

30° MICST™ 2.0 Angled Needle
BL3420A
RIS: 2.0mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



30° MicroFlow® Straight Needle
DP8230
RIS: 2.50-2.65mm OD: 1.07mm
ID: 0.90-0.51mm



45° MicroFlow® Plus Straight Needle
DP8345
RIS: 2.75-2.85mm OD:1.17mm
ID:0.98-0.61mm



30° MICST™ 2.0 Straight Needle
BL3420
RIS: 2.0mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



15° MicroFlow® Straight Needle
DP8215
RIS: 2.50-2.65mm OD: 1.07mm
ID: 0.90-0.51mm



30° MicroFlow® 2.2 Straight Needle
DP8730
RIS: 2.2mm OD: 0.90mm
ID: 0.70-0.50mm



30° MICST™ Angled Needle
BL3318A
RIS: 1.8-2.00mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



45° MicroFlow® Straight Needle
DP8245
RIS: 2.50-2.65mm OD: 1.07mm
ID: 0.90-0.51mm



15° Thin Tip 20G Straight Needle
DP8415
RIS: 2.5-2.65mm OD: 0.90mm
ID: 0.67mm



30° MICST™ Straight Needle
BL3318
RIS: 1.8-2.00mm OD: 0.74-0.95mm
ID: 0.5-0.79mm



30° MicroFlow® Plus Straight Needle
DP8330
RIS: 2.75-2.85mm OD:1.17mm
ID:0.98-0.61mm



30° Thin Tip 20G Straight Needle
DP8430
RIS: 2.5-2.65mm OD: 0.90mm
ID: 0.67mm



30° 19G Standard Straight Needle
DP8130
RIS: 3.00-3.20mm OD: 1.70mm
ID: 0.90mm



30° MicroFlow® Plus Angled Needle
DP8330A
RIS: 2.75-2.85mm OD:1.17mm
ID:0.98-0.61mm



30° Thin Tip 20G Angled Needle
DP8430A
RIS: 2.5-2.65mm OD: 0.90mm
ID: 0.67mm



Phaco Sleeves, Accessories and Kits

Kit with (1) BL3420S MICST™ 2.0mm Blue Straight Needle,
(2) DP5531 Green Sleeves, Wrench, Test Chamber and Plug
BLP324MX30
Sterile, 6 pouches/box



S/U Sleeve and Chamber Kit
1 Sleeve, 1 Chamber/Pouch
DP5501
Light blue standard D8109 sleeve and test chamber, sterile, 12 pouches/box



2.4mm Yellow Sleeve
BL3124
Limited use, non-sterile, 6/box



2.4mm Yellow Sleeve
2 Sleeves/Pouch
BL3124S
Yellow, sterile, 12 pouches/box



1.8mm MICST™ Sleeve
2 Sleeves/Pouch
BL3118
Translucent blue, sterile,
12 pouches /box



Thin Tip 20G Sleeve
2 Sleeves/Pouch
DP5511
Dark blue, sterile, 12 pouches/box



MicroFlow® 2.2mm Sleeve
2 Sleeves/Pouch
DP5531
Green, sterile, 12 pouches/box



Standard Sleeve
D8109
Light blue, limited use, non-sterile,
10/box





Irrigation/Aspiration

In order to complement the high levels of efficiency provided by Stellaris Elite®, we offer a range of irrigation/aspiration in both single use and reusable options. Our I/A handpieces, including our Capsuleguard and bimanual series, also provide the versatility to work with a variety of cataract platforms in addition to the dual-platform Stellaris Elite®.

With a large number of tips, materials and angles to choose from, surgeons from a wide range of practices with varied techniques can appreciate the numerous options available.



I/A
IRRIGATION/
ASPIRATION

I/A Single Use

Single-Use I/A MICST™
85795ST
I/A handpiece 45° angled tip, with sleeve for 1.8mm incision, 0.3mm aspiration port, single-use, 12/box



I/A Single Use CapsuleGuard®

CapsuleGuard® 45° I/A Handpiece
for all systems
85910S
Silicone tip, single-use, 12/box



CapsuleGuard® 45° I/A Handpiece for Stellaris®
85910ST
Silicone tip, single-use, 12/box



CapsuleGuard® Straight I/A Handpiece
for all systems
85912S
Silicone tip, single-use, 12/box



CapsuleGuard® Straight I/A Handpiece
for Stellaris®
85912ST
Silicone tip, single-use, 12/box



CapsuleGuard® Curved I/A Handpiece
for all systems
85913S
Silicone tip, single-use, 12/box



CapsuleGuard® Curved I/A Handpiece
for Stellaris®
85913ST
Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ Straight I/A Handpiece
for all systems
85914S
Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ Straight I/A Handpiece
for Stellaris®
85914ST
Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ 45° I/A Handpiece
for all systems
85915S
Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ 45° I/A Handpiece
for Stellaris®
85915ST
Silicone tip, single-use, 12/box



I/A Series BL3160 Handle

I/A Handle
BL3160
Stellaris® I/A handpiece with Surelock



Needle and sleeve not included.

Reusable I/A Tips for BL3160
BL3162
I/A tip 0.4mm straight silicone 20G



Reusable I/A Tips for BL3160
BL3163
I/A tip 0.4mm 45° angled, silicone 20G



Reusable I/A Tips for BL3160
BL3164
I/A tip 0.4mm straight metal 17.5G



Reusable I/A Tips for BL3160
BL3165
I/A tip 0.4mm curved metal 17.5G



Reusable I/A Tips for BL3160
BL3166
I/A tip 0.4mm 45° angled, metal 17.5G



I/A Series DP9721 Handle

I/A Handpiece
DP9721
I/A handpiece, titanium, luer connectors and knurled grip



Soft Tip I/A Straight
DP9739
Soft I/A tip 6/box Straight, limited reuse



Soft Tip I/A Curved
DP9740
Soft I/A tip 6/box Curved, limited reuse



Soft Tip I/A 45° Angled
DP9741
Soft I/A tip 6/box 45°, limited reuse



Reusable I/A Capsule Polisher Metal Tip
DP9744
I/A capsule polisher 45° Angled, Metal



Reusable I/A Capsule Polisher Tip Silicone Sleeve
DP9745
I/A capsule polisher 45° with silicone sleeve



Adaptor
DP9742
Adaptor for titanium handpiece



O-Rings
DP9728
Replacement rings for reusable I/A tips



I/A Tip Straight
DP9722
Fits handle DP9721, Straight I/A tip, 16G O.D. 0.3mm aspiration port



I/A Tip Curved
DP9723
Fits handle DP9721, Curved I/A tip, 16G O.D. 0.3mm aspiration port



I/A Tip 45° Angled
DP9724
Fits handle DP9721, Angled I/A tip, 16G 0.3mm aspiration port



I/A Tip Modified U Tip
DP9725
Fits handle DP9721, 180° Modified U tip, 16G, 0.3mm open end port



Aspirating Adaptor
DP9726
Male luer aspirating adaptor for use with aspirating cannulas, titanium




I/A Tip 90° Angled
DP9730
Fits handle DP9721, 90° Angled aspiration tip, 16G, 0.3mm aspiration port



I/A Tip 21G Straight

DP9731
Fits handle DP9721, Straight I/A tip, 21G, for use with silicone irrigation sleeve, 0.3mm aspiration port



I/A Tip 21G 45° Angled

DP9733
Fits handle DP9721, 45° Angled tip, 21G, for use with silicone irrigation sleeves, 0.3mm aspiration port



I/A Tip 21G Curved

DP9734
Fits handle DP9721, curved I/A tip, 21G, for use with silicone irrigation sleeve, 0.3mm aspiration port



I/A Handpiece Adapter 3" Tubing

MVS 1061L
Double lumen tubing, female connector on irrigation line, male on aspiration line



I/A Handpiece 17G Straight

MVS 1063
Straight 20mm shaft with tip guard, 0.3mm aspiration port, uses 75mm of double lumen tubing (MVS 1061L)



I/A Handpiece 17G 60°

MVS 1063C
Straight 20mm shaft, 2mm angled 60° tip with tip guard, 0.3mm aspiration port, uses 75mm of double lumen tubing (MVS 1061L)



I/A Handpiece 17G Curved

MVS 1063S
Curved 20mm shaft, with tip guard, 0.3mm aspiration port, uses 75mm of double lumen tubing (MVS 1061L)



Bimanual I/A Single Use - 21G Procedure for Stellaris® Systems

S/U Bimanual 21G Aspiration Handpiece

85780S
Roughened tip, 0.8mm O.D., Curved shaft, 12/box



S/U Bimanual 21G Irrigation Handpiece

85787ST
Smooth tip, 0.8mm O.D., Curved shaft, for Stellaris® only, 12/box



Bimanual I/A Single Use - 23G Procedure for Stellaris® Systems

S/U Bimanual 23G Aspiration Handpiece

85792S
Roughened tip, 0.6mm O.D., Curved shaft, 12/box



S/U Bimanual 23G Irrigation Handpiece

85807ST
Smooth tip, 0.6mm O.D. Curved shaft, for Stellaris® only, 12/box



S/U Bimanual 23G Irrigation Handpiece

85807S
Smooth tip, 0.6mm O.D. Curved shaft, 12/box



Bimanual I/A Single Use - 23G Procedure for other Phaco Systems

S/U Bimanual 23G Aspiration Handpiece
85792S
Roughened tip, 0.6mm O.D, Curved shaft, 12/box



S/U Bimanual 23G Irrigation Handpiece
85807S
SMOOTH tip, 0.6mm O.D, Curved shaft, 12/box



Reusable Bimanual

Bimanual 19G Irrigating Handpiece
MVS1082
Designed for small incision phaco, 0.6mm port, Titanium



Bimanual 18G Aspirating Handpiece
MVS1083
Designed for small incision phaco, .2mm roughened tip, 0.3mm port, Titanium



Bimanual 23G Aspiration Handpiece
ET4750 A23
0.4mm smooth tip, Titanium



Bimanual 23G Aspiration Handpiece
ET4750 AP23
0.4mm roughened tip, Titanium



Bimanual 21G Irrigating Handpiece
ET4750 I21
Dual 0.4mm irrigating ports, Titanium



Bimanual 23G Aspiration Handpiece
ET4753 A23
0.3mm smooth tip, Titanium



Instruments

With our large investment in R&D and manufacturing facilities, we not only have precision engineered single-use and reusable instruments for the most performed procedures, but also instruments that surgeons trust for rare and uncommon procedures. That is why so many have come to trust Bausch + Lomb and Storz® day in and day out, no matter what challenges arise.

Our dedication to evolving to meet the needs of the modern surgical practice is evident in our wide range of instruments, designed for reliable performance and trusted efficiency.



IN
INSTRUMENTS

Single-Use Knives

Single-Use 15° Stab Guarded Knife
LP7515G
15° sideport guarded, 6/box



Single-Use 1.5-2.0mm
Trapezoid Knife
LP7596
1.5-2.0mm trapezoidal, 6/box



Single-Use 15° Sideport Knife
LP7515
15° sideport, 6/box



Single-Use 2.75mm
Angled Slit Knife
LP7559A
2.75mm, angled up, 6/box



Single-Use 1.6-1.8mm
Trapezoid Knife
LP7600
1.6-1.8mm trapezoidal, 6/box



Single-Use 2.0-2.2mm
Trapezoid Knife
LP7601
2.0-2.2mm trapezoidal, 6/box



Single-Use 2.85mm
Angled Slit Knife
LP7553A
2.85mm, angled up, 6/box



Single-Use 1.8-2.2-2.8mm
Unitrapezoid Knife
LP7602
1.8-2.2-2.8mm uni-trap, 6/box



Single-Use 2.65mm
Angled Slit Knife
LP7551A
2.65mm, angled up, 6/box



Diamond Knives

Tru-Size™ Angled Trapezoid Diamond Knife
1.8-2.2mm
E0130 1.8/2.2



Diamond Profile Angled 2.65mm Keratome
E0130 2.65



Diamond Profile Angled 1.8mm Keratome
E0130 1.8MM



Diamond Profile Angled 2.7-2.9mm
Trapezoid Knife
E0130 2.75/2.9



Tru-Size™ Angled Trapezoid Diamond Knife
2.0-2.4mm
E0130 2.0/2.4



Diamond Profile Angled 2.8mm Keratome
E0130 2.8



Diamond Profile Angled 2.2mm Keratome
E0130 2.2MM



Reusable Capsulorhexis Forceps

Capsulorhexis Forceps Handle
ET8190 H
Capsulorhexis forceps handle titanium (non-rotating)



Capsulorhexis Forceps Tip 13mm Shaft
ET1089 T
For use with ET8190 H handle, for 1.4mm incision, curved shaft, 13mm shaft



Capsulorhexis Forceps Tip 15mm Shaft
ET1088 TS
For use with ET8190 H handle, for 1.4mm incision, curved shaft, 15mm shaft



Utrata MICST™ Capsulorhexis Forceps
E2002 M
Angled 13mm shafts, flat handle non-glare finish



Silverstein MICST™ 1.8mm Capsulorhexis Forceps
E2026
Ultra-thin profile, for sub-2.0mm incision, marks at 5.0mm and 6.0mm



Devine Micro Capsulorhexis Forceps
E2020 S
Angled sharp tips, curved shafts, round handle, for 1.8mm incision



Devgan Micro Capsulorhexis Forceps
ET2012
Ultra fine tip, cross-action forceps, marks at 2.5mm and 5.0mm



Utrata MICST™ Capsulorhexis Forceps
ET1202
13mm shaft, Titanium



Wiley Capsulorhexis Forceps
E2002 W
Angled 12mm shafts, iris stop, marks at 2.5mm and 5mm



Single-Use Capsulorhexis Forceps

Single-Use Capsulorhexis Forceps
SUE01-12
Straight jaws, fine sharp tips, 12/box



Single-Use MICST™ Capsulorhexis Forceps
SUE06-12
Fine angled tips for 1.8mm procedures, 12/box



Single-Use Marked Capsulorhexis Forceps
SUE07-12
Straight jaws with marks at 2.75-3.0mm and 5.5-6.0mm, 12/box



Single-Use MICST™ Curved Capsulorhexis Forceps
SUE08-12
MICST™ capsulorhexis forceps curved, 12/box



Accessories

With the dual-platform Stellaris Elite®, we have whole-eye solutions that are reliable and convenient. Surgeons trust Bausch + Lomb for their entire practice and a variety of procedures.

Our accessories are engineered to give reliable performance and are made with high-quality components. With easy reordering and robust customer service, we ensure that your practice is covered and make it easy for you to grow with Bausch + Lomb Surgical.



ACCESSORIES


Anterior Vitrectomy Packs - Stellaris®

Stellaris® Anterior Vitrectomy Pack
BL5612
20G vitrectomy cutter with male adaptor, tubing, infusion sleeve and irrigation cannula

- Need BL3190 Anterior Vitrectomy Adaptor Kit when using this pack with Stellaris®



20G Anterior Vitrectomy Cutter
MVS3201-6
20G vitrectomy cutter without irrigation



Anterior Vitrectomy - Stellaris® PC or Stellaris Elite® Anterior

Stellaris® Anterior Vitrectomy Pack
BL5612
20G vitrectomy cutter with male adaptor, tubing, infusion sleeve and irrigation cannula 6/box



23G Vitrectomy Cutter
BL5623S
23G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



20G Vitrectomy Cutter
BL5620
20G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



25G Vitrectomy Cutter
BL5625S
25G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



Anterior Vitrectomy - Stellaris Elite® Posterior

20G 7.5K Pneumatic Vitrectomy Cutter
BL5626
20G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



25G Stellaris Elite® Cutter
SE5625
25G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



23G Stellaris Elite® Cutter
SE5623
23G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



Anterior Vitrectomy

Anterior Vitrecomy Adapter Kit
BL3190
21G cannula, drive line, male adaptor, 6/box



Anterior Vitrectomy

Anterior Vitrecomy Adapter Kit

BL3190
21G cannula, drive line, male adaptor, 6/box



23G Single Valved Trocar

BL5283
For pars plana anterior vitrectomy



25G Single Valved Trocar

BL5285
For pars plana anterior vitrectomy



S/U Bimanual 22G Irrigation Handpiece

85806ST
Smooth tip, 0.7mm O.D., curved shaft, for Stellaris® only, 12/box



S/U Bimanual 23G Irrigation Handpiece

85807ST
Smooth tip, 0.6mm O.D., curved shaft, for Stellaris® only, 12/box



S/U Bimanual 21G Irrigation Handpiece


85787ST
Smooth tip, 0.8mm O.D., curved shaft for Stellaris® only, 12/box



Maintenance

Universal Maintenance Kit

CX7120
CX7123 Universal Aspiration Brush
CX7124 Universal 3-Way Stopcock
Kit contains syringe, stopcock and brush




Stellaris® PC Mayo Tray Adaptor

BL4364
Fits plastic arm for both Stellaris® PC and Stellaris ELITE®



DigiFlow™ Tubing

D4600A
Air exchange line



O-Rings

DP9728
Replacement rings for reusable I/A tips



Bipolar


Reusable Bipolar Cord
CX9430
with lemo connector



Reusable Bipolar Cord
S2050 B
2-pin connectors



Single-Use Bipolar Cord
S2070 C
2-pin connectors, 10/box



Single-Use Bipolar Pencil 18G Straight
E7918
2-pin connectors, 10/box



Single-Use Bipolar Pencil 18G Curved
E7918 C
2-pin connectors, 10/box



Banana Plug Adaptor
CX9404
reusable



Reusable Bipolar Cord
CX9400



Drapes/Covers

Front Panel Drape
DP5300
Clear plastic, 10/box

Mayo Stand Covers
DRA7352
Clear plastic, 24"X40", 50/box

Patient Half Body Incise Drape
DRA1061
Embossed plastic, one fluid pouch, sterile, 20/box

Patient Half-Body Aperture Drape
DRA1059A
Embossed plastic, one fluid pouch, sterile, 20/box

Patient 3/4 Length Aperture Drape
DRA830A
Blue SMS material, two fluid pouches, sterile, 20/box

Patient 3/4 Length Incise Drape
DRA831
Blue SMS material, two fluid pouches, sterile, 20/box



Standard Back Table Packs

Standard Cataract Procedure Pack
LYN023BSBLB
With full body aperture SMS drape, two fluid pouches and a separate incise film

Standard Vitrectomy Procedure Pack
LYN023VTBLC
With large EENT split drape and separate half body plastic incise drape

Basic Cataract Procedure Pack Incise
LYN08904SPC
With half body embossed plastic incise drape, one fluid pouch




Trays

Microsurgical Instrument Tray Small
E7409
Sterilizable plastic tray, capacity 4-7 instruments



Microsurgical Instrument Tray Double Layer
E7411
Sterilizable plastic tray, capacity 20-30 instruments



Microsurgical Instrument Tray Large
E7406
Sterilizable plastic tray, capacity 10-20 instruments



Stellaris® Ultrasonic Handpiece Tray
BL3379



Index

Stellaris Elite®		
BL2300HD	3	
BL3160	3	
BL3170	3	
D4600A	3	
DP9721	3	
Phaco Packs		
BL5111	6	
BL5112	7	
BL5113	7	
BL5113A	7	
BL5110	8	
BL5114	8	
BL5114A	8	
BLV2260000	9	
BLV3220000	9	
BLV4180000	9	
BL5115-2	10	
BL5115-4	10	
BL5115-1	11	
BL5115-3	11	
Phaco Needles and Sleeves		
BL3420AS	14	
BL3420S	14	
BL3318AS	14	
BL3318S	14	
DP8130S	14	
DP8145S	14	
DP8230S	14	
DP8245S	14	
DP8330S	14	
DP8430AS	14	
DP8430S	14	
DP8730S	14	
DP8930S	14	
BL3420A	15	
BL3420	15	
BL3318A	15	
BL3318	15	
DP8130	15	
DP8215	15	
DP8230	15	
DP8245	15	
DP8330	15	
DP8330A	15	
DP8345	15	
DP8415	15	
DP8430	15	
DP8430A	15	
DP8730	15	
BLP324MX30	16	
DP5501	16	
BL3124	17	
BL3124S	17	
BL3118	17	
DP5511	17	
DP5531	17	
D8109	17	
Irrigation/Aspiration		
85795ST	20	
85910S	20	
85910ST	20	
85912S	20	
85912ST	20	
85913S	21	
85913ST	21	
85914S	21	
85914ST	21	
85915S	21	
85915ST	21	
BL3160	22	
BL3162	22	
BL3163	22	
BL3164	22	
BL3165	22	
BL3166	22	
DP9728	23	
DP9721	23	
DP9739	23	
DP9740	23	
DP9741	23	
DP9744	23	
DP9745	23	
DP9742	23	
DP9722	23	
DP9723	23	
DP9724	23	
DP9725	23	
DP9726	23	
DP9730	23	
DP9731	24	
DP9733	24	
DP9734	24	
MVS 1061L	24	
MVS 1063	24	
MVS 1063C	24	
MVS 1063S	24	
85780S	25	
85787ST	25	
85792S	25	
85807ST	25	
85807S	25	
85792S	26	
85807S	26	
MVS1082	27	
MVS1083	27	
ET4750 A23	27	
ET4750 AP23	27	
ET4750 I21	27	
ET4753 A23	27	
Instruments		
LP7596	30	
LP7600	30	
LP7602	30	
LP7515G	30	
LP7515	30	
LP7601	30	
LP7551A	30	
LP7559A	30	
LP7553A	30	
E0130 1.8/2.2	31	
E0130 1.8mm	31	
E0130 2.0/2.4	31	
E0130 2.2mm	31	
E0130 2.65	31	
E0130 2.75/2.9	31	
E0130 2.8	31	
ET8190 H	32	
ET1089T	32	
ET1088TS	32	
E2002 M	32	
E2026	32	
E2020 S	32	
ET2012	32	
ET1202	32	
E2002 W	32	
SUE01-12	33	
SUE06-12	33	
SUE07-12	33	
SUE08-12	33	
Accessories		
BL5612	36	
MVS3201-6	36	
BL5612	36	
BL5620	36	
BL5623S	36	
BL5625S	36	
SE5623	37	
SE5625	37	
BL5626	37	
BL3190	38	
BL5283	38	
BL5285	38	
85806ST	38	
85807ST	38	
85787ST	38	
CX7120	39	
CX7123	39	
CX7124	39	
BL4364	39	
D4600A	39	
DP9728	39	
CX9430	40	
S2050 B	40	
S2070 C	40	
E7918	40	
E7918 C	40	
CX9404	40	
CX9400	40	
DP5300	41	
DRA7352	41	
DRA1061	41	
DRA1059A	41	
DRA830A	41	
DRA831	41	
LYN023BSBLB	42	
LYN023VTBLC	42	
LYN08904SPC	42	
E7409	43	
E7406	43	
E7411	43	
BL3379	43	

BAUSCH+LOMB SURGICAL SOLUTIONS

PARTNERING
IS OUR
PRACTICE.



For more information or to demo these technologies
contact your Bausch + Lomb representative at 800.338.2020

BAUSCH+LOMB



800.338.2020
stellariselite.com
customerservice@stellariselite.com

®/TM are trademarks of Bausch & Lomb Incorporated or its affiliates.
Any other product/brand names and/or logos are trademarks of the respective owners.
©2020 Bausch & Lomb Incorporated or its affiliates. SRG.0232.USA.19



WARNING: Class 4 Laser Product—visible laser radiation.
Avoid eye or skin exposure to direct or scattered radiation.